

ABX Minocal



LOT		CX 499		CAL		 (Exp.) 2024-11-05 (YYYY - MM - DD)					
Rev 1		ABX Lysebio									
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN	YUMIZEN								TOLERANCES TOLERANCE
		H1500	H1500								
		H2500	H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.40	8.40								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.50								± 0.06
HB HGB	g/dl	12.6	12.6								± 0.2
	g/l	126	126								± 2
	mmol/l	7.82	7.82								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	235								± 10
VGM MCV	µm ³ ; fl	79.1	79.1								± 4.0

1300116192
Ref. TEMP-0387 Rev.52