

ABX Minocal

LOT

CX 500

CAL



(Exp.) 2024-12-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9				8.7					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47				4.39					± 0.06
HB HGB	g/dl	13.1				13.0					± 0.2
	g/l	131				130					± 2
	mmol/l	8.14				8.07					± 0.12
HT HCT	%	37.3				37.3					± 1.0
	l/l	0.373				0.373					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	226				232					± 10
VMP MPV	μm ³ ; fl	8.7				9.5					± 0.5

ABX Minocal

LOT

CX 500

CAL



(Exp.) 2024-12-05

(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.5							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.40							± 0.06
HB HGB	g/dl		12.8							± 0.2
	g/l		128							± 2
	mmol/l		7.95							± 0.12
HT HCT	%		37.0							± 1.0
	l/l		0.370							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		229							± 10
VMP MPV *	µm ³ ; fl		8.7							± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300113555

* N/A for Micros Scil Vet ABC