

ABX Minocal



LOT		CX 500		CAL		(Exp.) 2024-12-05 (YYYY - MM - DD)						
		Rev 1		ABX Lysebio								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.9	9.0	8.9	8.9	8.7		8.9	8.9	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.46	4.39	4.46	4.45	4.44		4.48	4.46	± 0.06	
HB HGB	g/dl	12.9	13.0	13.0	13.2	12.9	13.0		13.2	13.2	± 0.2	
	g/l	129	130	130	132	129	130		132	132	± 2	
	mmol/l	8.01	8.07	8.07	8.20	8.01	8.07		8.20	8.20	± 0.12	
HT HCT	%	36.7	36.8	35.3	35.9	35.6	35.3		36.5	35.9	± 1.0	
	l/l	0.367	0.368	0.353	0.359	0.356	0.353		0.365	0.359	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	230	237	225	226	235	226		227	226	± 10	
VMP MPV	µm ³ ; fl	8.6	8.5	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

1300090308
BACK / VERSO
Ref: TEMP-0387 Rev.52