

ABX Minocal



LOT		CX 500		CAL		 (Exp.) 2024-12-05 (YYYY - MM - DD)					
Rev 1		ABX Lysebio									
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN	YUMIZEN								TOLERANCES TOLERANCE
		H1500	H1500								
		H2500	H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.35	8.35								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.52								± 0.06
HB HGB	g/dl	13.0	13.0								± 0.2
	g/l	130	130								± 2
	mmol/l	8.07	8.07								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	233	233								± 10
VGM MCV	µm ³ ; fl	82.2	82.2								± 4.0

1300116192
Ref. TEMP-0387 Rev.52