

ABX Minocal

LOT

CX 501

CAL



(Exp.) 2025-01-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9				8.8					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53				4.47					± 0.06
HB HGB	g/dl	13.2				13.1					± 0.2
	g/l	132				131					± 2
	mmol/l	8.20				8.14					± 0.12
HT HCT	%	37.4				37.8					± 1.0
	l/l	0.374				0.378					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	230				238					± 10
VMP MPV	μm ³ ; fl	8.5				9.1					± 0.5

ABX Minocal

LOT

CX 501

CAL



(Exp.) 2025-01-05

(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.5							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.45							± 0.06
HB HGB	g/dl		12.9							± 0.2
	g/l		129							± 2
	mmol/l		8.01							± 0.12
HT HCT	%		36.7							± 1.0
	l/l		0.367							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		243							± 10
VMP MPV *	µm ³ ; fl		8.3							± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300113555

* N/A for Micros Scil Vet ABC