

# ABX Minocal


 (Exp.) 2025-01-05  
 (YYYY - MM - DD)

| LOT                      |                                                       | CX 501      |         | CAL    |        | <div><div></div><div>(Exp.) 2025-01-05</div><div>(YYYY - MM - DD)</div></div> |        |        |                     |        |                         |
|--------------------------|-------------------------------------------------------|-------------|---------|--------|--------|-------------------------------------------------------------------------------|--------|--------|---------------------|--------|-------------------------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio |         |        |        |                                                                               |        |        |                     |        | TOLERANCES<br>TOLERANCE |
|                          |                                                       | MICROS      | MICROS  | PENTRA | PENTRA | PENTRA                                                                        | PENTRA | PENTRA | PENTRA              | PENTRA |                         |
|                          |                                                       | EMI CRP     | LC-767G | MS CRP | XLR    | DX NEXUS<br>DF NEXUS                                                          | MS60   |        | 60<br>60 C+<br>ES60 | XL80   |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.1         | 8.9     | 9.3    | 9.0    | 9.3                                                                           | 8.8    |        | 9.0                 | 9.0    | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.48        | 4.52    | 4.51   | 4.53   | 4.53                                                                          | 4.50   |        | 4.52                | 4.53   | ± 0.06                  |
| HB HGB                   | g/dl                                                  | 13.0        | 13.0    | 13.3   | 13.3   | 13.0                                                                          | 13.1   |        | 13.3                | 13.3   | ± 0.2                   |
|                          | g/l                                                   | 130         | 130     | 133    | 133    | 130                                                                           | 131    |        | 133                 | 133    | ± 2                     |
|                          | mmol/l                                                | 8.07        | 8.07    | 8.26   | 8.26   | 8.07                                                                          | 8.14   |        | 8.26                | 8.26   | ± 0.12                  |
| HT HCT                   | %                                                     | 36.7        | 36.8    | 36.3   | 36.2   | 35.6                                                                          | 35.6   |        | 36.6                | 36.2   | ± 1.0                   |
|                          | l/l                                                   | 0.367       | 0.368   | 0.363  | 0.362  | 0.356                                                                         | 0.356  |        | 0.366               | 0.362  | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 236         | 239     | 236    | 239    | 236                                                                           | 237    |        | 239                 | 239    | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 8.2         | 8.2     | N/A    | N/A    | N/A                                                                           | N/A    |        | N/A                 | N/A    | ± 0.5                   |

 1300090308  
 BACK / VERSO  
 Ref: TEMP-0387 Rev.52