

# ABX Minocal



**LOT**

**CX 501**

**CAL**



**(Exp.) 2025-01-05**  
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500	H1500	H2500	H2500					
		HELO 1.x	HELO 2.0							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.47</b>	<b>8.47</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.57</b>	<b>4.57</b>							± 0.06
HB HGB	g/dl	<b>13.2</b>	<b>13.2</b>							± 0.2
	g/l	<b>132</b>	<b>132</b>							± 2
	mmol/l	<b>8.20</b>	<b>8.20</b>							± 0.12
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>235</b>	<b>235</b>							± 10
VGM MCV	µm <sup>3</sup> ; fl	<b>80.7</b>	<b>80.7</b>							± 4.0

1300116192  
Ref: TEMP-0387 Rev.52