

ABX Minocal



LOT

CX 502

CAL



(Exp.) 2025-02-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9				8.8					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46				4.39					± 0.06
HB HGB	g/dl	13.2				13.2					± 0.2
	g/l	132				132					± 2
	mmol/l	8.20				8.20					± 0.12
HT HCT	%	37.7				37.3					± 1.0
	l/l	0.377				0.373					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	239				233					± 10
VMP MPV	µm ³ ; fl	8.0				7.1					± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300090308