

# ABX Minocal

**LOT**

**CX 502**

**CAL**



**(Exp.) 2025-02-05**  
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.6</b>				<b>8.7</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.39</b>				<b>4.41</b>					± 0.06
HB HGB	g/dl	<b>13.2</b>				<b>13.3</b>					± 0.2
	g/l	<b>132</b>				<b>133</b>					± 2
	mmol/l	<b>8.20</b>				<b>8.26</b>					± 0.12
HT HCT	%	<b>37.1</b>				<b>38.4</b>					± 1.0
	l/l	<b>0.371</b>				<b>0.384</b>					± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>224</b>				<b>235</b>					± 10
VMP MPV	μm <sup>3</sup> ; fl	<b>8.1</b>				<b>8.7</b>					± 0.5

# ABX Minocal

**LOT**

**CX 502**

**CAL**



(Exp.) **2025-02-05**  
(YYYY-MM-DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l		<b>8.1</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l		<b>4.32</b>							± 0.06
HB HGB	g/dl		<b>12.8</b>							± 0.2
	g/l		<b>128</b>							± 2
	mmol/l		<b>7.95</b>							± 0.12
HT HCT	%		<b>36.5</b>							± 1.0
	l/l		<b>0.365</b>							± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l		<b>235</b>							± 10
VMP MPV *	µm <sup>3</sup> ; fl		<b>8.0</b>							± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300113555

\* N/A for Micros Scil Vet ABC