

ABX Minocal



(Exp.) **2025-02-05**
(YYYY - MM - DD)

LOT

CX 502

CAL

Rev **1**

| PARAMETRES PARAMETERS | UNITES UNITS | WHITEDIFF | | | | | | | | | TOLERANCES TOLERANCE |
|--------------------------|---|----------------------------|--|--------------|-----------------|--------------|------------|---------------------|--------------------------------------|--------------|-------------------------|
| | | YUMIZEN | | YUMIZEN | YUMIZEN | YUMIZEN | | YUMIZEN | YUMIZEN | YUMIZEN | |
| | | H500 OT H500 CT H550 | | H500 OT | H500 CT H550 | H560 | | H500 OT H500E OT | H500 CT H500E CT H550 H550E | H560 | |
| v1.0 to v2.x | | v3.x Only | | | v3.0 to v3.3 | | Since v4.x | | | Since v3.4 | |
| GB WBC | 10 ³ /mm ³ ;10 ⁹ /l | 8.50 | | 8.50 | 8.50 | 8.50 | | 8.50 | 8.50 | 8.50 | ± 0.20 |
| GR RBC | 10 ⁶ /mm ³ ;10 ¹² /l | 4.43 | | 4.39 | 4.39 | 4.39 | | 4.41 | 4.41 | 4.41 | ± 0.06 |
| HB HGB | g/dl | 12.9 | | 12.9 | 12.9 | 12.9 | | 12.9 | 12.9 | 12.9 | ± 0.2 |
| | g/l | 129 | | 129 | 129 | 129 | | 129 | 129 | 129 | ± 2 |
| | mmol/l | 8.01 | | 8.01 | 8.01 | 8.01 | | 8.01 | 8.01 | 8.01 | ± 0.12 |
| HT HCT | % | 37.4 | | 36.5 | 35.7 | 35.7 | | 37.9 | 37.6 | 37.6 | ± 1.0 |
| | l/l | 0.374 | | 0.365 | 0.357 | 0.357 | | 0.379 | 0.376 | 0.376 | ± 0.010 |
| PLA PLT | 10 ³ /mm ³ ;10 ⁹ /l | 242 | | 235 | 246 | 246 | | 236 | 245 | 245 | ± 10 |
| VMP MPV | µm ³ ; fl | 9.6 | | N/A | N/A | N/A | | N/A | N/A | N/A | ± 0.5 |

Ref: TEMP-0387 Rev.52 BACK / VERSO 1300113555