

ABX Minocal


 (Exp.) 2025-02-05
 (YYYY - MM - DD)

LOT
CX 502
CAL

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN	YUMIZEN	YUMIZEN		YUMIZEN	YUMIZEN	YUMIZEN	
		H500 OT H500 CT H550		H500 OT	H500 CT H550	H560		H500 OT H500E OT	H500 CT H500E CT H550 H550E	H560	
		v1.0 to v2.x			v3.x Only	v3.0 to v3.3			Since v4.x	Since v3.4	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.50		8.50	8.50	8.50		8.50	8.50	8.50	± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.43		4.39	4.39	4.39		4.41	4.41	4.41	± 0.06
HB HGB	g/dl	12.9		12.9	12.9	12.9		12.9	12.9	12.9	± 0.2
	g/l	129		129	129	129		129	129	129	± 2
	mmol/l	8.01		8.01	8.01	8.01		8.01	8.01	8.01	± 0.12
HT HCT	%	37.4		36.5	35.7	35.7		37.9	37.6	37.6	± 1.0
	l/l	0.374		0.365	0.357	0.357		0.379	0.376	0.376	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	242		235	246	246		236	245	245	± 10
VMP MPV	μm ³ ; fl	9.6		N/A	N/A	N/A		N/A	N/A	N/A	± 0.5

 Ref: TEMP-0387 Rev.52
 BACK / VERSO
 1300113555