



Manufactured for
 HORIBA Instruments Incorporated:
 Pointe Brand
 5449 Research Drive
 Canton, MI 48188



LOT 421302

SAP 5390015524



REF T7535-CAL



2025-12-31

Calibrator Values*

	Valeur titrée Assigned value Deklariertes Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tildelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd	Unité Unit Einheit Unidad Unità Unidade Enhed Enhet Μονάδα Jednostka Enhhet
Calibrator A	0	mg/dL
	0.00	g/L
Calibrator B	82	mg/dL
	0.82	g/L
Calibrator C	164	mg/dL
	1.64	g/L
Calibrator D	273	mg/dL
	2.73	g/L
Calibrator E	382	mg/dL
	3.82	g/L
Calibrator F	546	mg/dL
	5.46	g/L

*Standardized against CAP/BCR/IFCC Reference Preparation for Proteins in Human Serum (RPPHS), lot 91/0619.