



Manufactured for
HORIBA Instruments Incorporated:
Pointe Brand
5449 Research Drive
Canton, MI 48188



LOT 401502

SAP 5390015523



REF F7505-CTL



2025-05-31

(FR) Constituant (EN) Constituent (DE) Bestandteil (ES) Componente (IT) Componente (PT) Componente (DA) Bestanddel (SV) Konstituent (EL) συστατικό (PL) Składnik (NO) Bestanddel	Unité Unit Einheit Unidad Unità Unidade Enhed Enhet Μονάδα Jednostka Enhet	Level I Exp. Date: 2025-05-31 Valeur titrée Assigned value Deklarierter Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tilldelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd			Level II Exp. Date: 2025-05-31 Valeur titrée Assigned value Deklarierter Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tilldelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd		
		Mean	SD	Range	Mean	SD	Range
Ferritin (1st I.S. Standardization)	ng/mL	105	8.0	89-121	410	30.8	349-472
Ferritin (3rd I.S. Standardization)	ng/mL	125	9.5	106-144	488	36.5	415-561

Ferritin values are provided with two options: 1st or 3rd I.S. (International Standard) Values.

It is recommended to use the same standardized system as was chosen when selecting values for the Pointe Ferritin Calibrators (F7505-CAL).