



**Yumizen G CTRL DDi I
 for YUMIZEN
 G1500/G1550**



D 4940630464140626D



D 5024074081D

Reference No.:		1300036414
Lot No.:		940630 (rev. 1)
Expiry date:		2026-06-30
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,24 – 0,74

**Yumizen G CTRL DDi I
 for YUMIZEN G800**



D 4940630464140626D



D 5024074081D

Reference No.:		1300036414
Lot No.:		940630 (rev. 1)
Expiry date:		2026-06-30
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,24 – 0,74

**Yumizen G CTRL DDi I
 for YUMIZEN
 G400 DDi/G200**

Reference No.:		1300036414
Lot No.:		940630 (rev. 1)
Expiry date:		2026-06-30
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,22 – 0,72

**Yumizen G CTRL DDi II
 for YUMIZEN
 G1500/G1550**



D 4 9 4 0 6 3 0 5 6 4 1 4 0 6 2 6 D



D 5 1 1 4 2 1 4 0 7 8 D

Reference No.:		1300036414
Lot No.:		940630 (rev. 1)
Expiry date:		2026-06-30
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,14 – 2,14

**Yumizen G CTRL DDi II
 for YUMIZEN G800**



D 4 9 4 0 6 3 0 5 6 4 1 4 0 6 2 6 D



D 5 1 1 3 2 1 3 0 7 6 D

Reference No.:		1300036414
Lot No.:		940630 (rev. 1)
Expiry date:		2026-06-30
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,13 – 2,13

**Yumizen G CTRL DDi II
 for YUMIZEN
 G400 DDi/G200**

Reference No.:		1300036414
Lot No.:		940630 (rev. 1)
Expiry date:		2026-06-30
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,30 – 2,30