



HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée
 B.P. 7290
 34184 Montpellier cedex 4
 FRANCE
 www.horiba.com



A03A00013C

Yumizen G FIB 5 for Yumizen G1500/G1550 and Yumizen G800



D 4 9 4 0 8 0 1 3 6 3 8 4 0 8 2 6 D



D 5 4 0 0 0 9 0 3 0 0 1 0 9 2 0 0 D



D 6 1 4 5 1 0 0 2 3 4 0 0 0 0 0 0 D



D 7 1 3 2 D

Reference No.:		1300036384	
Lot No.:		940801 (rev. 1)	
Expiry date:		2026-08-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550/G800	g/L	sec
		4	9,0
		3	10,9
		2	14,5
		1	23,4
Dilution		10x	

Yumizen G FIB 5 for Yumizen G400 DDi/G200



D 4 9 4 0 8 0 1 3 6 3 8 4 0 8 2 6 D



D 5 4 0 0 0 9 3 3 0 0 1 1 5 2 0 0 D



D 6 1 5 3 1 0 0 2 5 1 0 0 0 0 0 0 D



D 7 1 3 0 D

Reference No.:		1300036384	
Lot No.:		940801 (rev. 1)	
Expiry date:		2026-08-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G400 DDi/G200	g/L	sec
		4	9,3
		3	11,5
		2	15,3
		1	25,1
Dilution		10x	