

ABX Minocal

LOT

CX 503

CAL



(Exp.) 2025-03-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1				9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53				4.51					± 0.06
HB HGB	g/dl	13.3				13.4					± 0.2
	g/l	133				134					± 2
	mmol/l	8.26				8.32					± 0.12
HT HCT	%	37.4				38.6					± 1.0
	l/l	0.374				0.386					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	229				237					± 10
VMP MPV	μm ³ ; fl	8.4				9.2					± 0.5

ABX Minocal

LOT

CX 503

CAL



(Exp.) 2025-03-05

(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.50							± 0.06
HB HGB	g/dl		13.0							± 0.2
	g/l		130							± 2
	mmol/l		8.07							± 0.12
HT HCT	%		36.9							± 1.0
	l/l		0.369							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		233							± 10
VMP MPV *	µm ³ ; fl		8.3							± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300113555

* N/A for Micros Scil Vet ABC