

# ABX Minocal


 (Exp.) 2025-03-05  
 (YYYY - MM - DD)

| LOT                      |                                                       | CX 503  |         | CAL         |        | <div><div></div><div>(Exp.) 2025-03-05</div><div>(YYYY - MM - DD)</div></div> |        |        |                     |        |                         |  |
|--------------------------|-------------------------------------------------------|---------|---------|-------------|--------|-------------------------------------------------------------------------------|--------|--------|---------------------|--------|-------------------------|--|
|                          |                                                       | Rev 1   |         | ABX Lysebio |        |                                                                               |        |        |                     |        |                         |  |
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | MICROS  | MICROS  | PENTRA      | PENTRA | PENTRA                                                                        | PENTRA | PENTRA | PENTRA              | PENTRA | TOLERANCES<br>TOLERANCE |  |
|                          |                                                       | EMI CRP | LC-767G | MS CRP      | XLR    | DX NEXUS<br>DF NEXUS                                                          | MS60   |        | 60<br>60 C+<br>ES60 | XL80   |                         |  |
|                          |                                                       |         |         |             |        |                                                                               |        |        |                     |        |                         |  |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.4     | 9.3     | 9.6         | 9.2    | 9.4                                                                           | 9.1    |        | 9.3                 | 9.2    | ± 0.2                   |  |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.51    | 4.49    | 4.47        | 4.48   | 4.52                                                                          | 4.50   |        | 4.53                | 4.48   | ± 0.06                  |  |
| HB HGB                   | g/dl                                                  | 13.1    | 13.1    | 13.4        | 13.4   | 13.0                                                                          | 13.2   |        | 13.3                | 13.4   | ± 0.2                   |  |
|                          | g/l                                                   | 131     | 131     | 134         | 134    | 130                                                                           | 132    |        | 133                 | 134    | ± 2                     |  |
|                          | mmol/l                                                | 8.14    | 8.14    | 8.32        | 8.32   | 8.07                                                                          | 8.20   |        | 8.26                | 8.32   | ± 0.12                  |  |
| HT HCT                   | %                                                     | 36.8    | 36.8    | 36.2        | 36.3   | 35.9                                                                          | 36.0   |        | 36.7                | 36.3   | ± 1.0                   |  |
|                          | l/l                                                   | 0.368   | 0.368   | 0.362       | 0.363  | 0.359                                                                         | 0.360  |        | 0.367               | 0.363  | ± 0.010                 |  |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 235     | 239     | 220         | 231    | 240                                                                           | 231    |        | 235                 | 231    | ± 10                    |  |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | 8.2     | 8.5     | N/A         | N/A    | N/A                                                                           | N/A    |        | N/A                 | N/A    | ± 0.5                   |  |

 1300090308  
 BACK / VERSO  
 Ref: TEMP-0387 Rev.52