

ABX Minocal


 (Exp.) 2025-03-05
 (YYYY - MM - DD)

LOT		CX 503		CAL		(Exp.) 2025-03-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4	9.3	9.6	9.2	9.4	9.1		9.3	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.49	4.47	4.48	4.52	4.50		4.53	4.48	± 0.06
HB HGB	g/dl	13.1	13.1	13.4	13.4	13.0	13.2		13.3	13.4	± 0.2
	g/l	131	131	134	134	130	132		133	134	± 2
	mmol/l	8.14	8.14	8.32	8.32	8.07	8.20		8.26	8.32	± 0.12
HT HCT	%	36.8	36.8	36.2	36.3	35.9	36.0		36.7	36.3	± 1.0
	l/l	0.368	0.368	0.362	0.363	0.359	0.360		0.367	0.363	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	239	220	231	240	231		235	231	± 10
VMP MPV	µm ³ ; fl	8.2	8.5	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5

 1300090308
 BACK / VERSO
 Ref: TEMP-0387 Rev.52