



Yumizen G CTRL I (level 1) for Yumizen G Line



D 4940540464120526D



D 5107144078131079D



D 6118101137080134D



D 7077116083112101D



D 8168074110321435D



D 9291394192288229D



D 0343349D

Reference No.:		1300036412	
Lot No:		940540 (rev. 2)	
Expiry date:		2026-05-31	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	10,7-14,4
	Prothrombin %	%	78-131
	Prothrombin INR	INR	0,79-1,18
PT Liq	Prothrombin Time	sec	10,1-13,7
	Prothrombin %	%	80-134
	Prothrombin INR	INR	0,77-1,16
PT R	Prothrombin Time	sec	8,3-11,2
	Prothrombin %	%	101-168
	Prothrombin INR	INR	0,74-1,10
APTT	Partial Thromboplastin Time	sec	32,1-43,5
APTT Liq	Partial Thromboplastin Time	sec	29,1-39,4
FIB Cl	Fibrinogen	g/L	1,92-2,88
TT	Thrombin Time	sec	22,9-34,3
AT	Antithrombin	%	76-103

Important notice

The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.

The % and INR dimension are linked to this MNPT. Then, the ranges of the controls need to be adjusted.



Yumizen G CTRL II (level 2) for Yumizen G Line



D4940540564120526D



D5167251028047144D



D6216171256027044D



D7157235130196028D



D8047129193611826D



D9612828060141000D



D0999371D

Reference No.:		1300036412	
Lot No:		940540 (rev. 2)	
Expiry date:		2026-05-31	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	16,7-25,1
	Prothrombin %	%	28-47
	Prothrombin INR	INR	1,44-2,16
PT Liq	Prothrombin Time	sec	17,1-25,6
	Prothrombin %	%	27-44
	Prothrombin INR	INR	1,57-2,35
PT R	Prothrombin Time	sec	13,0-19,6
	Prothrombin %	%	28-47
	Prothrombin INR	INR	1,29-1,93
APTT	Partial Thromboplastin Time	sec	61,1-82,6
APTT Liq	Partial Thromboplastin Time	sec	61,2-82,8
FIB Cl	Fibrinogen	g/L	0,60-1,41
TT	Thrombin Time	sec	N/A
AT	Antithrombin	%	29-48

Important notice

The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.

The % and INR dimension are linked to this MNPT. Then, the ranges of the controls need to be adjusted.