



HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée
 BP 7290
 34184 Montpellier
 cedex 4
 FRANCE
 www.horiba.com



A03A00007B

Yumizen G FIB 2 for Yumizen G1500/G1550 and Yumizen G800



D 4 9 4 0 4 0 7 3 6 3 8 3 0 4 2 6 D



D 5 4 0 0 0 8 7 3 0 0 1 0 7 2 0 0 D



D 6 1 4 2 1 0 0 2 3 2 0 0 0 0 0 0 D



D 7 1 2 8 D

Reference No.:		1300036383	
Lot No.:		940407 (rev. 1)	
Expiry date:		2026-04-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550/G800	g/L	sec
		4	8,7
		3	10,7
		2	14,2
		1	23,2
	Dilution	10x	

Yumizen G FIB 2 for Yumizen G400 DDi/G200



D 4 9 4 0 4 0 7 3 6 3 8 3 0 4 2 6 D



D 5 4 0 0 0 9 4 3 0 0 1 1 3 2 0 0 D



D 6 1 4 5 1 0 0 2 2 3 0 0 0 0 0 0 D



D 7 1 2 6 D

Reference No.:		1300036383	
Lot No.:		940407 (rev. 1)	
Expiry date:		2026-04-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G400 DDi/G200	g/L	sec
		4	9,4
		3	11,3
		2	14,5
		1	22,3
	Dilution	10x	