



LAtrol™ Normal Control Plasma *LAtrol®*

HORIBA Reference No.:	1300079435
Biomedica Reference No.:	816N
Lot No.:	241001
Expiry date:	2026/10/01

for Yumizen G800

Parameter	Unit	Range
DRVVT	sec	43.7 – 56.1
DRVVC	sec	34.3 – 43.7
Normalized Ratio		<1.16

for Yumizen G1500/G1550

Parameter	Unit	Range
DRVVT	sec	32.7 – 56.1
DRVVC	sec	32.0 – 44.4
Normalized Ratio		<1.16

Important notice

Each laboratory should establish its own DRVVT and DRVVC normal reference range (mean value ± 2 S.D.).

The Normalized Ratio is calculated as follows:

$$\text{Normalized Ratio} = \frac{\text{DRVVT (sample)} / \text{DRVVT (mean clot time test)}}{\text{DRVVC (sample)} / \text{DRVVC (mean clot time confirm)}}$$