

ABX Minocal



LOT

CX 504

CAL



(Exp.) 2025-04-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0				8.9					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40				4.40					± 0.06
HB HGB	g/dl	13.0				13.0					± 0.2
	g/l	130				130					± 2
	mmol/l	8.07				8.07					± 0.12
HT HCT	%	37.0				37.0					± 1.0
	l/l	0.370				0.370					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	232				230					± 10
VMP MPV	µm ³ ; fl	8.3				7.4					± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300090308