

ABX Minocal


 (Exp.) 2025-04-05
 (YYYY - MM - DD)

LOT		CX 504		CAL		(Exp.) 2025-04-05 (YYYY - MM - DD)						
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE	
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA		
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	8.9	9.1	8.8	9.0	8.6		9.0	8.8	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.38	4.37	4.34	4.35	4.42	4.38		4.40	4.35	± 0.06	
HB HGB	g/dl	12.8	13.0	13.0	13.0	13.0	12.9		13.1	13.0	± 0.2	
	g/l	128	130	130	130	130	129		131	130	± 2	
	mmol/l	7.95	8.07	8.07	8.07	8.07	8.01		8.14	8.07	± 0.12	
HT HCT	%	36.4	36.5	35.2	35.5	35.6	35.0		36.1	35.5	± 1.0	
	l/l	0.364	0.365	0.352	0.355	0.356	0.350		0.361	0.355	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	230	232	225	226	231	228		232	226	± 10	
VMP MPV	µm ³ ; fl	8.3	8.3	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

 1300090308
 BACK / VERSO
 Ref: TEMP-0387 Rev.52