

# ABX Minocal



**LOT** CX 504  
Rev 1

**CAL**

(Exp.) 2025-04-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80	
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	9.0	8.9	9.1	8.8	9.0	8.6		9.0	8.8	± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	4.38	4.37	4.34	4.35	4.42	4.38		4.40	4.35	± 0.06
HB HGB	g/dl	12.8	13.0	13.0	13.0	13.0	12.9		13.1	13.0	± 0.2
	g/l	128	130	130	130	130	129		131	130	± 2
	mmol/l	7.95	8.07	8.07	8.07	8.07	8.01		8.14	8.07	± 0.12
HT HCT	%	36.4	36.5	35.2	35.5	35.6	35.0		36.1	35.5	± 1.0
	l/l	0.364	0.365	0.352	0.355	0.356	0.350		0.361	0.355	± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	230	232	225	226	231	228		232	226	± 10
VMP MPV	µm <sup>3</sup> ; fl	8.3	8.3	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.52 BACK / VERSO 1300090308