

ABX Minocal



LOT

CX 504

Rev 1

CAL



(Exp.) 2025-04-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500 HELO 1.x	H1500 H2500 HELO 2.0							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.30	8.30							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42							± 0.06
HB HGB	g/dl	13.1	13.1							± 0.2
	g/l	131	131							± 2
	mmol/l	8.14	8.14							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	235							± 10
VGM MCV	μm ³ ; fl	82.5	82.5							± 4.0

1300116192
Ref: TEMP-0387 Rev.52