

ABX Minocal



LOT		CX 504		CAL		 (Exp.) 2025-04-05 (YYYY - MM - DD)					
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H1500	H1500								
		H2500	H2500								
		HELO 1.x	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.30	8.30								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42								± 0.06
HB HGB	g/dl	13.1	13.1								± 0.2
	g/l	131	131								± 2
	mmol/l	8.14	8.14								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	235	235								± 10
VGM MCV	µm ³ ; fl	82.5	82.5								± 4.0

1300116192
Ref. TEMP-0387 Rev.52