

Yumizen G DDi 2 for YUMIZEN G1500/G1550



| | | | |
|----------------|---------------------|-----------------|--------|
| Reference No.: | | 1300036391 | |
| Lot No.: | | 940531 (rev. 1) | |
| Expiry date: | | 2025-05-31 | |
| Coagulometers | | MASTER CURVE | |
| Optical | Yumizen G1500/G1550 | µg(FEU)/mL | ΔOD |
| | | 4,00 | 0.2891 |
| | | 3,00 | 0.2220 |
| | | 2,00 | 0.1430 |
| | | 1,00 | 0.0650 |
| | | 0,50 | 0.0310 |
| | | 0,25 | 0.0150 |
| Slope value | | 0,07 | |

Yumizen G DDi 2 for YUMIZEN G800



| | | | |
|----------------|--------------|-----------------|--------|
| Reference No.: | | 1300036391 | |
| Lot No.: | | 940531 (rev. 1) | |
| Expiry date: | | 2025-05-31 | |
| Coagulometers | | MASTER CURVE | |
| Optical | Yumizen G800 | µg(FEU)/mL | ΔOD |
| | | 4,00 | 0.2720 |
| | | 3,00 | 0.2130 |
| | | 2,00 | 0.1390 |
| | | 1,00 | 0.0620 |
| | | 0,50 | 0.0290 |
| | | 0,25 | 0.0140 |
| Slope value | | 0,06 | |



Yumizen G DDi 2 for YUMIZEN G400 DDi/G200

| | | | |
|---|------------------------------|------------------------|---------------|
| Reference No.: | | 1300036391 | |
| Lot No.: | | 940531 (rev. 1) | |
| Expiry date: | | 2025-05-31 | |
| Coagulometers | | MASTER CURVE | |
| | | µg(FEU)/mL | ΔOD |
| Optical | Yumizen G400 DDi/G200 | 4,00 | 0.1140 |
| | | 3,00 | 0.0878 |
| | | 2,00 | 0.0595 |
| | | 1,00 | 0.0295 |
| | | 0,50 | 0.0145 |
| | | 0,25 | 0.0068 |
| | | Slope value | 0,02 |
| Important notice | | | |
| <p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type Lin-Lin p-p is selected in Calibration > D-Dimer > ug FEU/ml. For more information, please contact your local HORIBA Medical representative.</p> | | | |