

Yumizen G DDi 2 for YUMIZEN G1500/G1550



Reference No.:		1300036391	
Lot No.:		940915 (rev. 1)	
Expiry date:		2025-09-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550	µg(FEU)/mL	ΔOD
		4,00	0.2740
		3,00	0.2185
		2,00	0.1350
		1,00	0.0545
		0,50	0.0262
		0,25	0.0142
Slope value		0,07	

Yumizen G DDi 2 for YUMIZEN G800



Reference No.:		1300036391	
Lot No.:		940915 (rev. 1)	
Expiry date:		2025-09-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G800	µg(FEU)/mL	ΔOD
		4,00	0.2670
		3,00	0.2070
		2,00	0.1315
		1,00	0.0576
		0,50	0.0265
		0,25	0.0125
Slope value		0,06	



Yumizen G DDi 2 for YUMIZEN G400 DDi/G200

Reference No.:		1300036391	
Lot No.:		940915 (rev. 1)	
Expiry date:		2025-09-30	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4,00	0.1109
		3,00	0.0835
		2,00	0.0566
		1,00	0.0292
		0,50	0.0150
		0,25	0.0078
		Slope value	0,02
Important notice			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type Lin-Lin p-p is selected in Calibration > D-Dimer > ug FEU/ml. For more information, please contact your local HORIBA Medical representative.</p>			