

ABX Minocal



LOT

CX 505

CAL



(Exp.) 2025-05-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9				9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46				4.48					± 0.06
HB HGB	g/dl	13.1				13.2					± 0.2
	g/l	131				132					± 2
	mmol/l	8.14				8.20					± 0.12
HT HCT	%	37.2				37.4					± 1.0
	l/l	0.372				0.374					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	238				232					± 10
VMP MPV	µm ³ ; fl	7.7				7.1					± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300090308