

ABX Minocal


LOT
CX 505

 Rev **1**
CAL

(Exp.) 2025-05-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.8	9.2	9.0	9.2	8.9		8.9	9.0	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.46	4.43	4.43	4.47	4.50	4.45		4.48	4.47	± 0.06
HB HGB	g/dl	13.0	13.0	13.1	13.2	12.9	13.1		13.2	13.2	± 0.2
	g/l	130	130	131	132	129	131		132	132	± 2
	mmol/l	8.07	8.07	8.14	8.20	8.01	8.14		8.20	8.20	± 0.12
HT HCT	%	36.8	35.9	35.4	36.0	35.8	35.4		36.1	36.0	± 1.0
	l/l	0.368	0.359	0.354	0.360	0.358	0.354		0.361	0.360	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	232	233	227	227	232	231		230	227	± 10
VMP MPV	µm ³ ; fl	8.0	8.2	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.52 BACK / VERSO 1300090308