

**LOT**
**CX 505**
**CAL**

**(Exp.) 2025-05-05**  
 (YYYY - MM - DD)

 Rev **1**

| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX Lysebio  |              |              |              |                      |              |        |                     |              | TOLERANCES<br>TOLERANCE |
|--------------------------|-------------------------------------------------------|--------------|--------------|--------------|--------------|----------------------|--------------|--------|---------------------|--------------|-------------------------|
|                          |                                                       | MICROS       | MICROS       | PENTRA       | PENTRA       | PENTRA               | PENTRA       | PENTRA | PENTRA              | PENTRA       |                         |
|                          |                                                       | EMI CRP      | LC-767G      | MS CRP       | XLR          | DX NEXUS<br>DF NEXUS | MS60         |        | 60<br>60 C+<br>ES60 | XL80         |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | <b>9.1</b>   | <b>8.8</b>   | <b>9.2</b>   | <b>9.0</b>   | <b>9.2</b>           | <b>8.9</b>   |        | <b>8.9</b>          | <b>9.0</b>   | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | <b>4.46</b>  | <b>4.43</b>  | <b>4.43</b>  | <b>4.47</b>  | <b>4.50</b>          | <b>4.45</b>  |        | <b>4.48</b>         | <b>4.47</b>  | ± 0.06                  |
| HB HGB                   | g/dl                                                  | <b>13.0</b>  | <b>13.0</b>  | <b>13.1</b>  | <b>13.2</b>  | <b>12.9</b>          | <b>13.1</b>  |        | <b>13.2</b>         | <b>13.2</b>  | ± 0.2                   |
|                          | g/l                                                   | <b>130</b>   | <b>130</b>   | <b>131</b>   | <b>132</b>   | <b>129</b>           | <b>131</b>   |        | <b>132</b>          | <b>132</b>   | ± 2                     |
|                          | mmol/l                                                | <b>8.07</b>  | <b>8.07</b>  | <b>8.14</b>  | <b>8.20</b>  | <b>8.01</b>          | <b>8.14</b>  |        | <b>8.20</b>         | <b>8.20</b>  | ± 0.12                  |
| HT HCT                   | %                                                     | <b>36.8</b>  | <b>35.9</b>  | <b>35.4</b>  | <b>36.0</b>  | <b>35.8</b>          | <b>35.4</b>  |        | <b>36.1</b>         | <b>36.0</b>  | ± 1.0                   |
|                          | l/l                                                   | <b>0.368</b> | <b>0.359</b> | <b>0.354</b> | <b>0.360</b> | <b>0.358</b>         | <b>0.354</b> |        | <b>0.361</b>        | <b>0.360</b> | ± 0.010                 |
| PLA PLT                  | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | <b>232</b>   | <b>233</b>   | <b>227</b>   | <b>227</b>   | <b>232</b>           | <b>231</b>   |        | <b>230</b>          | <b>227</b>   | ± 10                    |
| VMP MPV                  | µm <sup>3</sup> ; fl                                  | <b>8.0</b>   | <b>8.2</b>   | N/A          | N/A          | N/A                  | N/A          |        | N/A                 | N/A          | ± 0.5                   |

 1300090308  
 BACK / VERSO  
 Ref: TEMP-0387 Rev.52