

# ABX Minocal



**LOT**

**CX 505**

Rev 1

**CAL**



**(Exp.) 2025-05-05**  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500	H1500	H2500	H2500					
		HELO 1.x	HELO 2.0							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.40</b>	<b>8.40</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.50</b>	<b>4.50</b>							± 0.06
HB HGB	g/dl	<b>13.1</b>	<b>13.1</b>							± 0.2
	g/l	<b>131</b>	<b>131</b>							± 2
	mmol/l	<b>8.14</b>	<b>8.14</b>							± 0.12
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>237</b>	<b>237</b>							± 10
VGM MCV	µm <sup>3</sup> ; fl	<b>81.0</b>	<b>81.0</b>							± 4.0

1300116192  
Ref: TEMP-0387 Rev.52