

Yumizen G FIB 2 for Yumizen G1500/G1550 and Yumizen G800



Reference No.:		1300036383	
Lot No.:		950113 (rev. 1)	
Expiry date:		2027-01-31	
Coagulometers		MASTER CURVE	
		g/L	sec
Optical	Yumizen G1500/G1550/G800	4	9,6
		3	11,7
		2	15,5
		1	25,2
		Dilution	10x

Yumizen G FIB 2 for Yumizen G400 DDi/G200



Reference No.:		1300036383	
Lot No.:		950113 (rev. 1)	
Expiry date:		2027-01-31	
Coagulometers		MASTER CURVE	
		g/L	sec
Optical	Yumizen G400 DDi/G200	4	9,8
		3	11,9
		2	15,6
		1	24,7
		Dilution	10x