

# ABX Minocal



**LOT**

**CX 506**

**CAL**



**(Exp.) 2025-06-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60 Care ST				60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.2</b>				<b>9.1</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.35</b>				<b>4.39</b>					± 0.06
HB HGB	g/dl	<b>13.1</b>				<b>13.1</b>					± 0.2
	g/l	<b>131</b>				<b>131</b>					± 2
	mmol/l	<b>8.14</b>				<b>8.14</b>					± 0.12
HT HCT	%	<b>37.2</b>				<b>37.1</b>					± 1.0
	l/l	<b>0.372</b>				<b>0.371</b>					± 0.010
PLA PLT	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>226</b>				<b>229</b>					± 10
VMP MPV	µm <sup>3</sup> ; fl	<b>8.0</b>				<b>7.2</b>					± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300090308