

ABX Minocal


 (Exp.) 2025-06-05
 (YYYY - MM - DD)

LOT		CX 506		CAL		(Exp.) 2025-06-05 (YYYY - MM - DD)						
		Rev 1		ABX Lysebio								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.0	9.3	9.1	9.3	9.2		9.3	9.1	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.37	4.37	4.34	4.36	4.40	4.34		4.39	4.36	± 0.06	
HB HGB	g/dl	12.9	13.0	13.0	13.2	13.0	12.9		13.2	13.2	± 0.2	
	g/l	129	130	130	132	130	129		132	132	± 2	
	mmol/l	8.01	8.07	8.07	8.20	8.07	8.01		8.20	8.20	± 0.12	
HT HCT	%	36.5	36.7	35.2	36.0	35.9	35.4		36.2	36.0	± 1.0	
	l/l	0.365	0.367	0.352	0.360	0.359	0.354		0.362	0.360	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	229	231	224	228	236	220		228	228	± 10	
VMP MPV	µm ³ ; fl	8.1	8.1	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

 1300090308
 BACK / VERSO
 Ref: TEMP-0387 Rev.52