

ABX Minocal

LOT CX 506		CAL		(Exp.) 2025-06-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H1500	H1500								
		H2500	H2500								
		HELO 1.X	HELO 2.0								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.60	8.60								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.43	4.43								± 0.06
HB HGB	g/dl	13.1	13.1								± 0.2
	g/l	131	131								± 2
	mmol/l	8.14	8.14								± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	234	234								± 10
VGM MCV	μm ³ ; fl	82.0	82.0								± 4.0

 1300116192
 Ref. TEMP-0387 Rev.52