

ABX Minocal

LOT
CX 506
CAL

(Exp.) 2025-06-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN	YUMIZEN	YUMIZEN		YUMIZEN	YUMIZEN	YUMIZEN	
		H500 OT H500 CT H550		H500 OT	H500 CT H550	H560		H500 OT H500E OT	H500 CT H500E CT H550 H550E	H560	
		v1.0 to v2.x			v3.x Only	v3.0 to v3.3			Since v4.x	Since v3.4	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85		8.85	8.85	8.85		8.85	8.85	8.85	± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.37		4.37	4.37	4.37		4.37	4.37	4.37	± 0.06
HB HGB	g/dl	12.7		12.8	12.8	12.8		12.7	12.7	12.7	± 0.2
	g/l	127		128	128	128		127	127	127	± 2
	mmol/l	7.89		7.95	7.95	7.95		7.89	7.89	7.89	± 0.12
HT HCT	%	36.5		35.0	34.5	34.5		36.5	37.0	37.0	± 1.0
	l/l	0.365		0.350	0.345	0.345		0.365	0.370	0.370	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	233		238	258	258		232	242	242	± 10
VMP MPV	μm ³ ; fl	9.8		N/A	N/A	N/A		N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.52 1300113555 BACK / VERSO