



Yumizen G CTRL DDi I for YUMIZEN G1500/G1550



D 4 9 4 0 9 2 1 4 6 4 1 4 0 9 2 6 D



D 5 0 2 8 0 7 8 0 9 5 D

Reference No.:		1300036414
Lot No.:		940921 (rev. 2)
Expiry date:		2026-09-30
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,28 – 0,78

Yumizen G CTRL DDi I for YUMIZEN G800



D 4 9 4 0 9 2 1 4 6 4 1 4 0 9 2 6 D



D 5 0 2 4 0 7 4 0 8 7 D

Reference No.:		1300036414
Lot No.:		940921 (rev. 2)
Expiry date:		2026-09-30
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,24 – 0,74

Yumizen G CTRL DDi I for YUMIZEN G400 DDi/G200

Reference No.:		1300036414
Lot No.:		940921 (rev. 2)
Expiry date:		2026-09-30
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,21 – 0,71



Yumizen G CTRL DDi II for YUMIZEN G1500/G1550



D 4 9 4 0 9 2 1 5 6 4 1 4 0 9 2 6 D



D 5 1 0 4 2 0 4 0 8 2 D

Reference No.:		1300036414
Lot No.:		940921 (rev. 2)
Expiry date:		2026-09-30
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,04 – 2,04

Yumizen G CTRL DDi II for YUMIZEN G800



D 4 9 4 0 9 2 1 5 6 4 1 4 0 9 2 6 D



D 5 0 9 7 1 9 7 1 0 4 D

Reference No.:		1300036414
Lot No.:		940921 (rev. 2)
Expiry date:		2026-09-30
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,97 – 1,97

Yumizen G CTRL DDi II for YUMIZEN G400 DDi/G200

Reference No.:		1300036414
Lot No.:		940921 (rev. 2)
Expiry date:		2026-09-30
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,03 – 2,03