

ABX Minocal



LOT

CX 507

CAL



(Exp.) 2025-07-05

(YYYY - MM - DD)

Rev 1

| PARAMETRES PARAMETERS | UNITES UNITS | ABX Lyse | | | | | | | | TOLERANCES TOLERANCE | |
|--------------------------|---|-----------------|--------|--------|--------|-----------------|--------|--------|--------|-------------------------|---------|
| | | MICROS | MICROS | MICROS | MICROS | MICROS | PENTRA | PENTRA | PENTRA | | PENTRA |
| | | ES60 Care ST | | | | 60 8P 60 18P | | | | | |
| GB WBC | 10 ³ /mm ³ ;10 ⁹ /l | 9.1 | | | | 9.1 | | | | | ± 0.2 |
| GR RBC | 10 ⁶ /mm ³ ;10 ¹² /l | 4.45 | | | | 4.45 | | | | | ± 0.06 |
| HB HGB | g/dl | 13.2 | | | | 13.3 | | | | | ± 0.2 |
| | g/l | 132 | | | | 133 | | | | | ± 2 |
| | mmol/l | 8.20 | | | | 8.26 | | | | | ± 0.12 |
| HT HCT | % | 37.4 | | | | 37.4 | | | | | ± 1.0 |
| | l/l | 0.374 | | | | 0.374 | | | | | ± 0.010 |
| PLA PLT | 10 ³ /mm ³ ;10 ⁹ /l | 231 | | | | 234 | | | | | ± 10 |
| VMP MPV | µm ³ ; fl | 7.9 | | | | 7.2 | | | | | ± 0.5 |

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300090308