

ABX Minocal


 (Exp.) 2025-07-05
 (YYYY - MM - DD)

LOT		CX 507		CAL		(Exp.) 2025-07-05 (YYYY - MM - DD)						
		Rev 1		ABX Lysebio								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.1	9.2	9.1	9.3	9.0		9.0	9.1	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.43	4.46	4.40	4.42	4.49	4.42		4.45	4.42	± 0.06	
HB HGB	g/dl	13.1	13.1	13.3	13.2	13.1	13.1		13.2	13.2	± 0.2	
	g/l	131	131	133	132	131	131		132	132	± 2	
	mmol/l	8.14	8.14	8.26	8.20	8.14	8.14		8.20	8.20	± 0.12	
HT HCT	%	36.8	36.8	36.5	36.5	35.9	35.8		36.5	36.5	± 1.0	
	l/l	0.368	0.368	0.365	0.365	0.359	0.358		0.365	0.365	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	234	239	231	229	238	229		236	229	± 10	
VMP MPV	µm ³ ; fl	8.1	8.2	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

 1300090308
 BACK / VERSO
 Ref: TEMP-0387 Rev.52