

ABX Minocal


LOT
CX 507
CAL

(Exp.) 2025-07-05
 (YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.1	9.2	9.1	9.3	9.0		9.0	9.1	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.43	4.46	4.40	4.42	4.49	4.42		4.45	4.42	± 0.06
HB HGB	g/dl	13.1	13.1	13.3	13.2	13.1	13.1		13.2	13.2	± 0.2
	g/l	131	131	133	132	131	131		132	132	± 2
	mmol/l	8.14	8.14	8.26	8.20	8.14	8.14		8.20	8.20	± 0.12
HT HCT	%	36.8	36.8	36.5	36.5	35.9	35.8		36.5	36.5	± 1.0
	l/l	0.368	0.368	0.365	0.365	0.359	0.358		0.365	0.365	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	234	239	231	229	238	229		236	229	± 10
VMP MPV	µm ³ ; fl	8.1	8.2	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5

 1300090308
 BACK / VERSO
 Ref: TEMP-0387 Rev.52