



Yumizen G PT Liq 4 for Yumizen G Line



Reference No.:		1300036373		
Lot No.:		940908 (rev. 1)		
Expiry date:		2026-03-31		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		1,22	%
		100		12,2
		50		17,2
		33		22,4
		25		27,2
	MNPT	12,2		
Important notice				
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p>				