



Yumizen G PT Liq 4 for Yumizen G Line



| | | | | |
|---|-----------------------|------------------------|---------------------|-------------|
| Reference No.: | | 1300036373 | | |
| Lot No.: | | 941115 (rev. 1) | | |
| Expiry date: | | 2026-05-31 | | |
| Coagulometers | | ISI | MASTER CURVE | |
| Optical | Yumizen G Line | | 1,26 | % |
| | | 100 | | 11,7 |
| | | 50 | | 17,2 |
| | | 33 | | 22,8 |
| | | 25 | | 28,0 |
| | | | MNPT | 11,7 |
| Important notice | | | | |
| The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT. | | | | |