

ABX Minocal

LOT

CX 508

CAL



(Exp.) 2025-08-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60				60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.7				9.2					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.59				4.59					± 0.06
HB HGB	g/dl	13.7				13.9					± 0.2
	g/l	137				139					± 2
	mmol/l	8.51				8.63					± 0.12
HT HCT	%	38.6				39.2					± 1.0
	l/l	0.386				0.392					± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	222				230					± 10
VMP MPV	μm ³ ; fl	8.0				8.3					± 0.5

ABX Minocal

LOT

CX 508

CAL



(Exp.) **2025-08-05**

(YYYY-MM-DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS	MICROS							
			Scil Vet ABC Plus Scil Vet ABC ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		8.7							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.54							± 0.06
HB HGB	g/dl		13.6							± 0.2
	g/l		136							± 2
	mmol/l		8.45							± 0.12
HT HCT	%		38.6							± 1.0
	l/l		0.386							± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l		223							± 10
VMP MPV *	µm ³ ; fl		7.9							± 0.5

Ref: TEMP-0387 Rev.52 FRONT / RECTO 1300113555

* N/A for Micros Scil Vet ABC