

ABX Minocal


 (Exp.) 2025-08-05
 (YYYY - MM - DD)

LOT		CX 508		CAL		(Exp.) 2025-08-05 (YYYY - MM - DD)							
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lysebio							TOLERANCES TOLERANCE		
				MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA			PENTRA
				EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60				60 60 C+ ES60
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.0	9.3	9.0	9.3	8.9		8.9	9.0	± 0.2	
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.55	4.50	4.55	4.61	4.53		4.61	4.55	± 0.06	
HB	HGB	g/dl	13.6	13.6	13.7	13.8	13.7	13.6		13.8	13.8	± 0.2	
		g/l	136	136	137	138	137	136		138	138	± 2	
		mmol/l	8.45	8.45	8.51	8.57	8.51	8.45		8.57	8.57	± 0.12	
HT	HCT	%	38.2	38.0	36.9	37.5	37.1	36.7		37.8	37.5	± 1.0	
		l/l	0.382	0.380	0.369	0.375	0.371	0.367		0.378	0.375	± 0.010	
PLA	PLT	10 ³ /mm ³ ;10 ⁹ /l	230	233	227	229	231	225		230	229	± 10	
VMP	MPV	µm ³ ; fl	7.9	8.1	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

 1300090308
 BACK / VERSO
 Ref: TEMP-0387 Rev.52