

Yumizen G DDi 2 for Yumizen G15X0x and Yumizen G405



Reference No.:		1300036391	
Lot No.:		950513 (rev. 1)	
Expiry date:		2026-05-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G15X0x/G405	µg(FEU)/mL	ΔOD
		4.00	0.2740
		3.00	0.2010
		2.00	0.1240
		1.00	0.0570
		0.50	0.0280
		0.25	0.0148
Slope value		0,07	

Yumizen G DDi 2 for Yumizen G8X0x



Reference No.:		1300036391	
Lot No.:		950513 (rev. 1)	
Expiry date:		2026-05-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G8X0x	µg(FEU)/mL	ΔOD
		4.00	0.2675
		3.00	0.2100
		2.00	0.1295
		1.00	0.0540
		0.50	0.0255
		0.25	0.0130
Slope value		0,06	



Yumizen G DDi 2 for Yumizen G400 DDi/G200

Reference No.:		1300036391	
Lot No.:		950513 (rev. 1)	
Expiry date:		2026-05-31	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4.00	0.1152
		3.00	0.0885
		2.00	0.0560
		1.00	0.0255
		0.50	0.0130
		0.25	0.0073
		Slope value	0,02
Important notice			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type Lin-Lin p-p is selected in Calibration > D-Dimer > ug FEU/ml. For more information, please contact your local HORIBA Medical representative.</p>			