

ABX Minocal


 (Exp.) 2025-09-05
 (YYYY - MM - DD)

LOT
CX 509
CAL

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2	9.4	9.0	9.3	9.0		9.0	9.0	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.55	4.50	4.54	4.59	4.54		4.59	4.54	± 0.06
HB HGB	g/dl	13.4	13.2	13.3	13.4	13.3	13.2		13.5	13.4	± 0.2
	g/l	134	132	133	134	133	132		135	134	± 2
	mmol/l	8.32	8.20	8.26	8.32	8.26	8.20		8.38	8.32	± 0.12
HT HCT	%	37.6	37.3	36.0	36.8	36.5	36.1		37.2	36.8	± 1.0
	l/l	0.376	0.373	0.360	0.368	0.365	0.361		0.372	0.368	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	230	232	233	229	240	230		235	229	± 10
VMP MPV	µm ³ ; fl	7.9	8.0	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5

 1300090308
 BACK / VERSO
 Ref: TEMP-0387 Rev.52