

ABX Minocal


 (Exp.) 2025-09-05
 (YYYY - MM - DD)

LOT		CX 509		CAL		(Exp.) 2025-09-05 (YYYY - MM - DD)						
		Rev 1		ABX Lysebio								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
		EMI CRP	LC-767G	MS CRP	XLR	DX NEXUS DF NEXUS	MS60		60 60 C+ ES60	XL80		
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2	9.4	9.0	9.3	9.0		9.0	9.0	± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.55	4.50	4.54	4.59	4.54		4.59	4.54	± 0.06	
HB HGB	g/dl	13.4	13.2	13.3	13.4	13.3	13.2		13.5	13.4	± 0.2	
	g/l	134	132	133	134	133	132		135	134	± 2	
	mmol/l	8.32	8.20	8.26	8.32	8.26	8.20		8.38	8.32	± 0.12	
HT HCT	%	37.6	37.3	36.0	36.8	36.5	36.1		37.2	36.8	± 1.0	
	l/l	0.376	0.373	0.360	0.368	0.365	0.361		0.372	0.368	± 0.010	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	230	232	233	229	240	230		235	229	± 10	
VMP MPV	µm ³ ; fl	7.9	8.0	N/A	N/A	N/A	N/A		N/A	N/A	± 0.5	

 1300090308
 BACK / VERSO
 Ref: TEMP-0387 Rev.52