

ABX Minocal

LOT		CX 509		CAL		(Exp.) 2025-09-05 (YYYY - MM - DD)						
Rev 1		ABX Lysebio										
PARAMETRES PARAMETERS	UNITES UNITS	YUMIZEN	YUMIZEN								TOLERANCES TOLERANCE	
		H1500	H1500									
		H2500	H2500									
		HELO 1.x	HELO 2.0									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.55	8.55								± 0.20	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.56								± 0.06	
HB HGB	g/dl	13.4	13.4								± 0.2	
	g/l	134	134								± 2	
	mmol/l	8.32	8.32								± 0.12	
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	232	232								± 10	
VGM MCV	μm ³ ; fl	82.0	82.0								± 4.0	

1300116192

Ref: TEMP-0387 Rev.52

 1300116192
 Ref. TEMP-0387 Rev.52