

LOT
CX 509
CAL

(Exp.) 2025-09-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN	YUMIZEN	YUMIZEN		YUMIZEN	YUMIZEN	YUMIZEN	
		H500 OT H500 CT H550 v1.0 to v2.x		H500 OT v3.x Only	H500 CT H550	H560 v3.0 to v3.3		H500 OT H500E OT Since v4.x	H500 CT H500E CT H550 H550E	H560 Since v3.4	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85		8.70	8.70	8.70		8.70	8.70	8.70	± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.63		4.51	4.51	4.51		4.56	4.56	4.56	± 0.06
HB HGB	g/dl	13.0		12.9	12.9	12.9		12.9	12.9	12.9	± 0.2
	g/l	130		129	129	129		129	129	129	± 2
	mmol/l	8.07		8.01	8.01	8.01		8.01	8.01	8.01	± 0.12
HT HCT	%	37.7		35.6	34.7	34.7		37.6	37.3	37.3	± 1.0
	l/l	0.377		0.356	0.347	0.347		0.376	0.373	0.373	± 0.010
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	240		249	256	256		242	255	255	± 10
VMP MPV	μm ³ ; fl	9.9		N/A	N/A	N/A		N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.52 1300113555 BACK / VERSO